

MS15X-S2

Precision XY Linear Stage



Overview

Multi-functional X/Y movement and rotational base designed for controlled adjustment in inspection providing fine movement with X and Y axis. The precision machine crafted components provide for smooth movements and adjustments in positions.

The MS15X-S2 can be conveniently attached to compatible stands (RK-10A, RK-10FS, RK-06A, RK-06FA, MS35B, MS37B, MS45BX-D2) with the use of C-Clamps or the more secure option: the MS15X-P4 mounting bracket that is included in the packaging. Also included is a Rotating Table that can screw directly on top of the MS15X base, altogether this allows for full range of movement in rotational, X and Y movements for the absolute best and precise inspection with its fine adjustment knobs.

Travel Area: X:1.5" / Y:1.25"

X/Y Table Size: 3.25" x 4.25"

MS15X-S1 travels ~18mm per full turn

MS15X-S2(R2) travels ~10mm per full turn

Note: Current version update MS15X-S2(R2) – August 2024

Specification

Model: MS15X-S2

Unit Weight: 907g (1.99lb)

Unit Dimension: 12.7cm (W) x 15.2cm (L) x 7.6cm (H)

<https://www.dinolite.us>

Toll-free: 888-668-2442

Local: 310-618-9990

International: +1 310 618 9990

Sales: sales@dinolite.us

Support: support@dinolite.us

Schedule a free online demo!

<https://www.dinolite.us/meeting>