

RK-10FS

High-precision Stand with Antistatic Baseplate



Overview

The RK-10FS is an upgraded version of the RK-10, which incorporates a heavier antistatic baseplate that suppresses external vibrations, making it an optimal choice for stable operation including high-magnification inspections. The baseplate rubber pads help make this an ESD dissipative stand.


The updated RK-10FS features an improved quick release button with an anti-backlash mechanism that enhances its ability to make precise focus adjustments, ensuring smoother and more responsive operation.

Specification

Model: RK-10FS

Unit Weight: 1.9kg

Unit Dimension: 22.6 cm (H) x 15 cm (D) x 22 cm (L)



<https://www.dinolite.us>

Toll-free: 888-668-2442

Local: 310-618-9990

International: +1 310 618 9990

Sales: sales@dinolite.us

Support: support@dinolite.us

Schedule a free online demo!

<https://www.dinolite.us/meeting>