

RK-10

Precision Stand with Quick Release



Overview

The Dino-Lite RK-10 universal stand is a premium, highly stable tabletop platform designed to give you complete command over your microscope's orientation. Fully compatible with the entire line of Dino-Lite handheld microscope cameras, it maximizes both productivity and ease of use.

- **Precision Focusing:** A dedicated fine-focus adjustment knob provides the exact control needed for sharp images, which is especially crucial when working at high magnifications.
- **Rapid Positioning:** A convenient quick-release button allows for fast, effortless vertical movement to accommodate different subjects.
- **Durable Construction:** Expertly crafted from resilient stainless steel and lightweight aluminum for long-lasting stability without unnecessary bulk.
- **Expandable Versatility:** Easily customize your workstation with optional attachments. The [RK-10-VX](#) can be added for extended vertical working distance, integrate seamlessly with the [RK-10-EX](#) horizontal arm for advanced reach and rotation, or the [RK-10-PX2](#) precision arm for smooth, micro-controlled XY target positioning.

PRODUCT SPECIFICATIONS

Model Number	RK-10		
Dimensions	22 cm (		
Weight	1050g (		
Body Material	Stand: Metal, Holster: Plastic		
Compatible with	All Dino-Lite Handheld Microscopes		
Package Includes	Stand, Mounting Arm, Holster		
Warranty Period	1 year		
Manufactured in	Taiwan		
Regulatory Approval	ROHS		

<https://www.dinolite.us>

Toll-free: 888-668-2442

Local: 310-618-9990

International: +1 310 618 9990

Sales: sales@dinolite.us

Support: support@dinolite.us

Schedule a free online demo!

<https://www.dinolite.us/meeting>